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ITL4
CE 2008
EN 149:2001 + A1:2009
REGULATION (UE) 2016/425

www.realcaresrl.it

(EN) RISK THAT THE PPE IS DESIGNED TO PROTECT FROM:

The KARINA ITL4 model is a FFP2 NR filtering half-mask. It is therefore designed to protect the user's respiratory tract against a potentially polluted atmosphere. The device is designed to provide protection against both solid and liquid aerosols.

(EN) CHECK BEFORE USE:

Before use, inspect the device for holes, tears or damage. If there is any damage, DO NOT USE. Get a new one.

(EN) HOW TO WEAR THE HALF MASK:

1. Hold the half mask so that the rounded outer part rests on the palm of the hand and that the elastics remain under the hand.
2. Place the half mask under the chin with the underwire laced up.
3. Place the rubber bands behind the ears.
4. Model the underwire nose clip to the shape of the nose itself. Check the seal of the half mask before entering the work area with the seal test, inhaling and exhaling deeply. If you feel air leaks, reposition the half mask.

(EN) STORAGE:

Store and transport unused half masks inside their packaging in a dry place (<80%) at a temperature between -20°C and +40°C. Disposable half mask: the dust filter half mask must not be used for more than one work shift (8 hours); do not reuse. Do not use after the expiration date indicated on the package. The head harness has the same manufacturing and expiry date as the complete device. Do not dispose of the device in the environment after use. Dispose of in accordance with local waste management regulations.

(EN) WARNING:

Check the integrity of the half mask before use. Replace the half mask immediately if you experience difficulty in breathing or if the half mask is damaged, or if there are problems with adherence to the face. The beard can reduce the efficiency of the half mask. Do not use when contaminant concentrations are hazardous to health or life or in oxygen deficiency. Do not use in an atmosphere where there is a danger of explosion. It does not require maintenance.

(EN) REGULATORY REFERENCES:


The CE marking certifies that the product meets the essential requirements of Reg. (EU) 2016/425 on personal protective equipment (PPE).

(EN) DECLARATION OF CONFORMITY:

Declaration of conformity downloadable from www.realcaresrl.it/wp-content/uploads/2021/04/FT01-DPI-ITL4-dichiarazione-di-conformita_Rev.2.pdf

(EN) PERFORMANCE RECORDED DURING TECHNICAL TESTS:

Total inward leakage:
less than 11% for at least 46 out of 50 tests performed; less than 8% for at least 8 out of 10 arithmetic averages for each wearer
Penetration of the filtering material:
Sodium chloride aerosol test 95 l/min <6%
Paraffin oil aerosol test 95 l/min <6%
Carbon dioxide content of the inhalation air: < 1.0%.
Respiratory resistance
Inhalation 30 l/min <0.7mbar
Inhalation 95 l/min <2.4mbar
Exhalation 160 l/min <3.0mbar
The test results qualify the device as a FFP2 NR class anti-dust half-mask.

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